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## Die Autoren

Rodolphe Barrangou  
Ph. D. Fred Breidt  
NC State University  
Department of Food Science  
322 Schaub Hall, Box 7624  
USA-Raleigh, NC 27695

Ph. D. Daniel Schroen  
Kendro Laboratory Products  
Applications Department  
31 Pecks Lane  
USA-Newtown CT 06470.

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# BANDING OF INTACT BACTERIOPHAGE USING THE S100-AT6 ROTOR AND DISCOVERY™ M150 MICRO-ULTRACENTRIFUGE

Rodolphe Barrangou and Fred Breidt, Ph.D., Department of Food Science,  
322 Schaub Hall, Box 7624, NC State University, Raleigh, NC 27695

Daniel Schroen, Ph.D., Applications Department, Kendro Laboratory Products,  
31 Pecks Lane, Newtown, CT 06470.

SORVALL® Discovery M150 Micro-ultracentrifuge

## Abstract

With recent advances in molecular biology, the isolation of bacteriophage continues to increase. The SORVALL Discovery M150 micro-ultracentrifuge allows rapid concentration and purification of bacteriophage in a CsCl step gradient. The pure phage solution obtained via ultracentrifugation can be subsequently used for electron microscopy and protein or nucleic acid analysis.

## Introduction

Phage are bacterial viruses which are ubiquitous in nature, and are likely to be found in the natural environment of their host organism.

They have been extensively studied in the food fermentation industry since they were identified as a main source of fermentation failure, especially in applications based on starter cultures (Jarvis, 1989; Hill, 1993; Brussow, 1998).

They have also been isolated from several natural environments, including soil (Cresswell, 1992; Pantastico-

Caldas, 1992; Ashelford, 1999; Ashelford et al., 2000), water, lakes and oceans (Bergh et al., 1989; Bratbak et al., 1990; Bratbak et al., 1994; Hennes and Simon, 1995; Wichels, 1998), and plants. Phage also interact with bacteria naturally present in those environments.

In order to characterize a particular bacteriophage, it must first be isolated and purified. It is therefore necessary to develop and optimize protocols that enable a rapid and efficient phage extraction.

Centrifugation has long been considered a convenient and relevant method to purify and concentrate bacteriophage. In this brief, we describe a rapid and convenient centrifugation protocol which can be used to purify phage. The phage concentrate obtained can be used directly for electron microscopy and nucleic acid purification.

## Materials and methods

A 100 ml phage lysate was prepared prior to phage purification by challenging 100 mL of a mid-log *Pediococcus* culture grown in MRS medium containing 10 mM CaCl<sub>2</sub> and incubated until clearing. The phage supernatant was collected by Superspeed centrifugation (SORVALL GSA or equivalent rotor, 20 min., 8,000 rpm, 10,400  $\times g$ ), and filter-sterilized using a 0.45 micron filtration membrane (Corning). The phage lysate was incubated overnight with 2.9% NaCl and 10% PEG at 4°C to promote phage precipitation. Phage were subsequently pelleted by centrifugation (8,000 rpm, 20 min., SORVALL GSA rotor), and the pellet was resuspended in 2 mL of TE buffer (100 mM Tris, pH 7.6, 50 mM EDTA, Sigma).

The 2 mL concentrated phage suspension was overlaid onto a three-step CsCl gradient containing 1 mL of 1.7 g/mL CsCl (molecular biology grade, Sigma), 1 mL of 1.5 g/mL CsCl and 1 mL of 1.4 g/mL CsCl in a 5 mL SORVALL ultracentrifuge tube.

Phage were centrifuged for 7 hr at 600,000  $\times g$  using the S100-AT6 rotor in a Discovery™ M150 Micro-ultracentrifuge. Phage-containing bands (translucent white/gray) were

extracted through the wall of the centrifuge tube by puncturing with a needle, and the CsCl was subsequently removed by dialysis (6,000 - 8,000 dalton membrane, Baxter Diagnostics Inc., McGraw Park, IL) for 15 hr with three changes of de-ionized water. Phage purification methods were mostly adapted from Maniatis et al. 1982. The resulting phage samples can be further used for electron microscopy or nucleic acid purification.

For electron microscopy, the concentrated phage were negatively stained with 2% uranyl acetate (pH 4.0) and electron microscopy was carried out using a JEOL model 100S electron microscope (JEOL USA, Peabody, MA), at 80 KV.

## Results

Phage purification and concentration using the CsCl step centrifugation yielded highly concentrated samples (over 10<sup>11</sup> phage / ml) of pure bacteriophage. The main advantage of this centrifugation step is time efficiency. In this case, only 7 hours of centrifugation with a Micro-ultracentrifuge were necessary, whereas 24 hours of centrifugation are typically required with a standard ultracentrifuge. Also, the samples obtained were highly concentrated, and very clean, enabling direct use for electron microscopy and protein or nucleic acid purification. Electron micrographs obtained with *Pediococcus* phage isolated from an industrial vegetable fermentation are shown in Figure 1.


Figure 1: Electron micrographs of *Pediococcus* bacteriophages. The phage particles are negatively stained with uranyl acetate at a magnification of 85,000 x. A, *oPs05a*; B, *oPs05b*; C, *oPs08a*; D, *Ps08b*; E, *Ps10a*; F, *Ps10b*. Each bar denotes 100 nm.

## For Ordering or Technical Information

Asia Pacific North	Kendro Laboratory Products (H.K.) Limited · Hong Kong · Tel +852 2711-3910 · Fax +852 2711-3858 · <a href="mailto:info@kendro.com">info@kendro.com</a>
Asia Pacific South	Kendro Laboratory Products · Sydney · NSW 2066 · Tel +61 (2) 9936-1540 · Fax +61 (2) 9427 9765 · <a href="mailto:info@kendro.com.au">info@kendro.com.au</a>
Europe, Middle East, Africa	Kendro Laboratory Products International Sales · Hanau · Germany · Tel +49 (1805) 536 370 · Fax +49 (1805) 112 114 · <a href="mailto:info@kendro.de">info@kendro.de</a>
Latin America	Kendro Laboratory Products International Sales · Hanau · Germany · Tel +49 (1805) 536 370 · Fax +49 (1805) 112 114 · <a href="mailto:info@kendro.de">info@kendro.de</a>
USA	Overseas Contract Services, Inc. · Wilmington, DE · USA · Tel +1 (302) 658-9272 · Fax +1 (302) 658-9970 · <a href="mailto:ocsKaren@delanet.com">ocsKaren@delanet.com</a>
Internet	Kendro Laboratory Products · Newtown · Tel +1 (800) 522-7746 · Fax +1 (203) 270-2166 · <a href="mailto:info@kendro.com">info@kendro.com</a> <a href="http://www.kendro.com">http://www.kendro.com</a>

Quality Products – Lifetime Care

Registered to ISO 9001. Kendro Laboratory Products meet or exceed stringent quality and product safety standards: CE for the European Union, and UL, cUL, or CSA standards for North America. ©2001 Kendro Laboratory Products. All Rights Reserved.

## Conclusions

The centrifugation procedure described in this study is designed for rapid and convenient phage purification and concentration. The suggested protocol for this application is 7 hours at 600,000 x g using the S100-AT6 rotor in a Discovery M150 Micro-Ultracentrifuge. The phage samples prepared using this method can be used for electron microscopy, and protein or nucleic acid purification.

## References

- Ashelford, K. E., M. J. Day, M. J. Bailey, A. K. Lilley, and J. C. Fry. 1999. In situ population dynamics of bacterial viruses in a terrestrial environment. *Appl. Environ. Microbiol.* 65:169-174.
- Ashelford, K. E., S. J. Norris, J. C. Fry, M. J. Bailey and M. J. Day. 2000. Seasonal population dynamics and interactions of competing bacteriophages and their host in the rhizosphere. *Appl. Environ. Microbiol.* 66:4193-4199.
- Bergh, O., Y. Borsheim, G. Bratbak, and M. Heldal. 1989. High abundance of viruses found in aquatic environments. *Nature*. 340:467-468.
- Bratbak, G., M. Heldal, S. Norland, and T. F. Thingstad. 1990. Viruses as partners in spring bloom microbial trophodynamics. *Appl. Environ. Microbiol.* 56:1400-1405.
- Bratbak, G., F. Thingstad, and M. Heldal. 1994. Viruses and the microbial loop. *Microb. Ecol.* 28:209-221.
- Brusow, H., A. Bruttin, F. Desiere, S. Lucchini, and S. Foley. 1998. Molecular ecology and evolution of *Streptococcus thermophilus* bacteriophages – a review. *Virus Genes*. 16:95-109.
- Cresswell, N., P. R. Herron, V. A. Saunders, and E. M. H. Wellington. 1992. The fate of introduced streptomycetes, plasmid and phage populations in a dynamic soil system. *J. Gen. Microbiol.* 138:659-666.
- Hennes, K. P., and M. Simon. 1995. Significance of bacteriophages for controlling bacterioplankton growth in a mesotrophic lake. *Appl. Environ. Microbiol.* 61:333-340.
- Hill, C. 1993. Bacteriophage and bacteriophage resistance in lactic acid bacteria. *FEMS Microbiol. Rev.* 12:87-108.
- Jarvis, A. W. 1989. Bacteriophages of lactic acid bacteria. *J. Dairy. Sci.* 72:3406-3428.
- Maniatis, T., E. F. Fritsch, and J. Sambrook. 1982. *Molecular Cloning: a laboratory manual*. Cold Spring Harbor Laboratory, Cold Spring Harbor, N. Y.
- Pantastico-Caldas, M., K. E. Duncan, C. A. Istock, and J. A. Bell. 1992. Population dynamics of bacteriophage and *Bacillus subtilis* in soil. *Ecology*. 73:1888-1902.
- Wichels, A., S. F. Biel, H. R. Gelderblom, T. Brinkhoff, G. Muyzer, and C. Schutt. 1998. Bacteriophage diversity in the North Sea. *Appl. Environ. Microbiol.* 64:4128-4133.

SORVALL®