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### Die Autoren

Rodolphe Barrangou  
Ph. D. Fred Breidt  
NC State University  
Department of Food Science  
322 Schaub Hall, Box 7624  
USA-Raleigh, NC 27695

Ph. D. Daniel Schroen  
Kendro Laboratory Products  
Applications Department  
31 Pecks Lane  
USA-Newtown CT 06470.

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Easy Info Nr. • 687

extracted through the wall of the centrifuge tube by puncturing with a needle; and the CsCl was subsequently removed by dialysis (6,000 - 8,000 dalton membrane, Baxter Diagnostics Inc., McGaw Park, IL) for 15 hr with three changes of de-ionized water. Phage purification methods were mostly adapted from Maniatis et al. 1982. The resulting phage samples can be further used for electron microscopy or nucleic acid purification.

For electron microscopy, the concentrated phage were negatively stained with 2% uranyl acetate (pH 4.0) and electron microscopy was carried out using a JEOL model 100S electron microscope (JEOL USA, Peabody, MA), at 80 KV.

## Results

Phage purification and concentration using the CsCl step centrifugation yielded highly concentrated samples (over  $10^{11}$  phage / ml) of pure bacteriophage. The main advantage of this centrifugation step is time efficiency. In this case, only 7 hours of centrifugation with a Micro-ultracentrifuge were necessary, whereas 24 hours of centrifugation are typically required with a standard ultracentrifuge. Also, the samples obtained were highly concentrated, and very clean, enabling direct use for electron microscopy and protein or nucleic acid purification. Electron micrographs obtained with *Pediococcus* phage isolated from an industrial vegetable fermentation are shown in Figure 1.


Figure 1: Electron micrographs of *Pediococcus* bacteriophages. The phage particles are negatively stained with uranyl acetate at a magnification of 85,000 x. A,  $\phi$ Ps05a; B,  $\phi$ Ps05b; C,  $\phi$ Ps08a; D, Ps08b; E, Ps10a; F, Ps10b. Each bar denotes 100 nm.

## Conclusions

The centrifugation procedure described in this study is designed for rapid and convenient phage purification and concentration. The suggested protocol for this application is 7 hours at 600,000 x g using the S100-AT6 rotor in a Discovery M150 Micro-Ultracentrifuge. The phage samples prepared using this method can be used for electron microscopy, and protein or nucleic acid purification.

## References

- Ashelford, K. E., M. J. Day, M. J. Bailey, A. K. Lilley, and J. C. Fry. 1999. In situ population dynamics of bacterial viruses in a terrestrial environment. *Appl. Environ. Microbiol.* 65:169-174.
- Ashelford, K. E., S. J. Norris, J. C. Fry, M. J. Bailey and M. J. Day. 2000. Seasonal population dynamics and interactions of competing bacteriophages and their host in the rhizosphere. *Appl. Environ. Microbiol.* 66:4193-4199.
- Bergh, O., Y. Borsheim, G. Bratbak, and M. Heidal. 1989. High abundance of viruses found in aquatic environments. *Nature.* 340:467-468.
- Bratbak, G., M. Heidal, S. Norland, and T. F. Thingstad. 1990. Viruses as partners in spring bloom microbial trophodynamics. *Appl. Environ. Microbiol.* 56:1400-1405.
- Bratbak, G., F. Thingstad, and M. Heidal. 1994. Viruses and the microbial loop. *Microb. Ecol.* 28:209-221.
- Brussow, H., A. Bruttin, F. Desiere, S. Lucchini, and S. Foley. 1998. Molecular ecology and evolution of *Streptococcus thermophilus* bacteriophages - a review. *Virus Genes.* 16:95-109.
- Cresswell, N., P. R. Herron, V. A. Saunders, and E. M. H. Wellington. 1992. The fate of introduced streptomycetes, plasmid and phage populations in a dynamic soil system. *J. Gen. Microbiol.* 138:659-666.
- Hennes, K. P., and M. Simon. 1995. Significance of bacteriophages for controlling bacterioplankton growth in a mesotrophic lake. *Appl. Environ. Microbiol.* 61:333-340.
- Hill, C. 1993. Bacteriophage and bacteriophage resistance in lactic acid bacteria. *FEMS Microbiol. Rev.* 12:87-108.
- Jarvis, A. W. 1989. Bacteriophages of lactic acid bacteria. *J. Dairy. Sci.* 72:3406-3428.
- Maniatis, T., E. F. Fritsch, and J. Sambrook. 1982. *Molecular Cloning: a laboratory manual.* Cold Spring Harbor Laboratory, Cold Spring Harbor, N. Y.
- Pantastico-Caldas, M., K. E. Duncan, C. A. Istock, and J. A. Bell. 1992. Population dynamics of bacteriophage and *Bacillus subtilis* in soil. *Ecology.* 73:1888-1902.
- Wichels, A., S. F. Biel, H. R. Gelderblom, T. Brinkhoff, G. Muyzer, and C. Schutt. 1998. Bacteriophage diversity in the North Sea. *Appl. Environ. Microbiol.* 64:4128-4133.

## For Ordering or Technical Information

<b>Asia Pacific North</b>	Kendro Laboratory Products (H.K.) Limited - Hong Kong - Tel: +852 2711-3970 - Fax: +852 2711-3858 - info@kendro.com
<b>Asia Pacific South</b>	Kendro Laboratory Products - Sydney NSW 2066 - Tel: +61 (2) 9936 1540 - Fax: +61 (2) 9427 9765 - info@kendro.com.au
<b>Europe, Middle East, Africa</b>	Kendro Laboratory Products International Sales - Hanau - Germany - Tel: +49 (1805) 536 376 - Fax: +49 (1805) 112 114 - info@kendro.de
<b>Latin America</b>	Overseas Contract Services, Inc. - Wilmington, DE, USA - Tel: +1 (302) 658-9272 - Fax: +1 (302) 658-9970 - ocs@kendro.com
<b>USA</b>	Kendro Laboratory Products - Newtown - Tel: +1 (800) 522-7746 - Fax: +1 (203) 270-2166 - info@kendro.com
<b>Internet</b>	http://www.kendro.com

Quality Products - Lifetime Care

SORVALL