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Die Autoren

Rodolphe Barrangou
Ph. D. Fred Breidt
NC State University
Department of Food Science
322 Schaub Hall, Box 7624
USA-Raleigh, NC 27695

Ph. D. Daniel Schroen
Kendro Laboratory Products
Applications Department
31 Pecks Lane
USA-Newtown CT 06470.

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Easy Info Nr. • 687

BANDING OF INTACT BACTERIOPHAGE USING THE S100-AT6 ROTOR AND DISCOVERY™ M150 MICRO-ULTRACENTRIFUGE

Rodolphe Barrangou and Fred Breidt, Ph.D., Department of Food Science,
322 Schaub Hall, Box 7624, NC State University, Raleigh, NC 27695

Daniel Schroen, Ph.D., Applications Department, Kendro Laboratory Products,
31 Pecks Lane, Newtown, CT 06470.

SORVALL® Discovery M150 Micro-ultracentrifuge

Abstract

With recent advances in molecular biology, the isolation of bacteriophage continues to increase. The SORVALL Discovery M150 micro-ultracentrifuge allows rapid concentration and purification of bacteriophage in a CsCl step gradient. The pure phage solution obtained via ultracentrifugation can be subsequently used for electron microscopy and protein or nucleic acid analysis.

Introduction

Phage are bacterial viruses which are ubiquitous in nature, and are likely to be found in the natural environment of their host organism.

They have been extensively studied in the food fermentation industry since they were identified as a main source of fermentation failure, especially in applications based on starter cultures (Jarvis, 1989; Hill, 1993, Brussow, 1998).

They have also been isolated from several natural environments, including soil (Cresswell, 1992; Pantastico-

Caldas, 1992; Ashelford, 1999; Ashelford et al., 2000), water, lakes and oceans (Bergh et al., 1989; Bratbak et al., 1990; Bratbak et al., 1994; Hennes and Simon, 1995; Wichels, 1998), and plants. Phage also interact with bacteria naturally present in those environments.

In order to characterize a particular bacteriophage, it must first be isolated and purified. It is therefore necessary to develop and optimize protocols that enable a rapid and efficient phage extraction.

Centrifugation has long been considered a convenient and relevant method to purify and concentrate bacteriophage. In this brief, we describe a rapid and convenient centrifugation protocol which can be used to purify phage. The phage concentrate obtained can be used directly for electron microscopy and nucleic acid purification.

Materials and methods

A 100 ml phage lysate was prepared prior to phage purification by challenging 100 mL of a mid-log *Pediococcus* culture grown in MRS medium containing 10 mM CaCl₂ and incubated until clearing. The phage supernatant was collected by Superspeed centrifugation (SORVALL GSA or equivalent rotor, 20 min., 8,000 rpm, 10,400 x g), and filter-sterilized using a 0.45 micron filtration membrane (Coming). The phage lysate was incubated overnight with 2.9% NaCl and 10% PEG at 4°C to promote phage precipitation. Phage were subsequently pelleted by centrifugation (8,000 rpm, 20 min., SORVALL GSA rotor), and the pellet was resuspended in 2 mL of TE buffer (100 mM Tris, pH 7.6, 50 mM EDTA, Sigma).

The 2 mL concentrated phage suspension was overlaid onto a three-step CsCl gradient containing 1 mL of 1.7 g/mL CsCl (molecular biology grade, Sigma), 1 mL of 1.5 g/mL CsCl and 1 mL of 1.4 g/mL CsCl in a 5 mL SORVALL ultracentrifuge tube.

Phage were centrifuged for 7 hr at 600,000 x g using the S100-AT6 rotor in a Discovery™ M150 Micro-ultracentrifuge. Phage-containing bands (translucent white/gray) were