

STAINING FOR PERCENTAGE AM COLONIZATION

([INVAM](#), [Giovannetti and Mosse, 1980](#))

Materials

250 or 500 um sieve

10% KOH

2% HCl

bunsen burner

beakers

tissue cassettes, vials or test tubes

0.05% Direct Blue stain [33 ml lactic acid, 33 ml glycerol, 33 ml water, 0.5 g Direct Blue 15 (Sigma D2535)]

Methods

- 1) Wash roots with tap water over sieve to remove soil.
- 2) Remove a subsample to determine dry weight.
- 3) Remove another subsample and place in tissue cassettes, vials or test tubes for staining.
- 4) Roots may be stored in water at 4°C like this for a week or so before staining. (If they need to be stored for a longer period, place the roots in a 1:1:1 lactic acid: glycerol: water solution at 4°C.)
- 5) Bring 10% KOH solution to a boil and pour over roots. Incubate for 5 min.
- 6) Remove KOH and rinse with tap water 5 times.
- 7) Add RT 2% HCl and incubate for 5 min.
- 8) Meanwhile, bring 0.05% Direct Blue stain to boil.
- 9) Remove HCl solution and add stain. Incubate 5 min.
- 10) Remove stain and store in water at 4°C. (Note: Stain is a hazardous waste, dispose of properly.)
- 11) If a long term storage is desired, store in 1:1:1 acid:glycerol:water solution at 4°C in air tight containers. Monitor for fungal contamination and if roots become destained, repeat steps 7-10.
- 12) Determine percentage colonization via grid-line intersect method outlined in [Giovannetti and Mosse, 1980](#).