The background of the cover is a close-up photograph of a metallic, cylindrical component, likely a sample introduction system for an ICP-MS, with a glass tube at the top. The lighting creates a strong rainbow-like iridescent effect across the metallic surface.

Spectroscopy[®]

Solutions for Materials Analysis

February 2008 Volume 23 Number 2

www.spectroscopyonline.com

Pittcon 2008

ICP-MS for Pharmaceutical Analysis

SERS Basics

Miniature Optical Spectrometers, Part II

Label-Free SERS for Rapid Species Identification of *Escherichia Coli*, *Listeria Monocytogenes*, and *Salmonella Typhimurium* Bacteria

The need for rapid detection of foodborne bacteria has long been a hot topic for a wide range of individuals, from policy makers to manufacturers. The traditional method, nucleic acid-based polymerase chain reaction (PCR), and antibody-based biosensors have been developed as viable tools to identify the bacteria. Generally, these methods are labor-intensive and multistep procedures, and they also require species-specific antibodies–probe reagents. Colloidal silver surface-enhanced Raman scattering (SERS), as demonstrated here, could be an important alternative technique in the rapid and simultaneous screening of the presence of the three most common outbreak bacteria because the technique uses exclusive biomarkers, is label-free, and provides easy sampling.

Yongliang Liu, Yud-Ren Chen, Xiangwu Nou, Moon S. Kim, and Kuanglin Chao



Frequent recalls of bacteria-contaminated ready-to-eat food products are great concerns for federal regulators and manufacturers, not only because of the public health issue but also because of the subsequent economic loss. To meet the requirement for rapid detection of foodborne bacteria, a number of new techniques, such as nucleic acid-based polymerase chain reaction (PCR) and capture antibody–detection (label) antibody-based sandwich immunoassays (1,2), have been developed. However, these methods have several potential problems, such as a high rate of false negatives and false positives, slow speed, selection of species-specific antibodies and probe reagents, and multistep sample preprocessing with the use of numerous chemical reagents.

Fast microbial detection requires minimal sample preparation, permits routine analysis with a shorter data collec-

tion interval and simpler operation, and it is hoped, reduces the probability of bacterial cross-contamination during the identification process. Surface-enhanced Raman scattering (SERS) studies of bacteria adsorbed on fresh borohydride-reduced or citrate-reduced silver colloids is an alternative approach. Not only has it been used to obtain the fingerprinting characterization of bacterial structure (3–5), but it also has an easy sampling attribute by mixing silver colloids directly with incubated cultures (5). We have identified the unique SERS bands in *Escherichia coli* and *Listeria monocytogenes* bacteria and subsequently have developed simple two-band algorithms for rapid and routine identification of *E. coli* and *L. monocytogenes* cultures on silver colloidal nanoparticles that were prepared from different batches and had undergone various storage durations (5).

Step up

with the proven MB3000
and its new features and applications



HORIZON MB™
FTIR Software

Ease of use and versatility.

Get the most bang for your FTIR buck

The most reliable FTIR system.
The lowest cost of ownership.


Learn more by visiting ABB's team
in booth # 6451

MB3000

Always ready. Every day.

ABB
Analytical Business Unit
Email: ftir@ca.abb.com
www.abb.com/analytical
www.mb3000ftir.com

The Right Choice For Optical Fiber



Fiberguide Industries produces up to 50,000 meters of optical fiber each day. Our four fiber drawing towers enable us to complete large orders on time and our extensive product line assures that we will have the right fiber for your application. Fiberguide's specialty optical fibers range from our extreme multi-mode silica/silica fiber Superguide G — which is suitable for space exploration or jet aircraft — to our graded index fiber suitable for a LAN.

Fiberguide's optical fibers are made with a variety of jackets depending on their use. With our specialty drawing tower we can apply metal jackets. Our extrusion capabilities allow us to apply any number of protective polymeric jackets to the fiber. Fiberguide's entire manufacturing process has strict quality control guidelines based on ISO 9000:2000 standards.

Fiberguide's high quality product and our exemplary customer service make us the premier manufacturer of specialty optical fiber. Fiberguide — the right choice for optical fiber.

Contact Fiberguide today and ask for our new CD-ROM. It's free!



1 Bay Street
Stirling, NJ 07980
Telephone: 908-647-6601
Fax: 908-647-8464
info@fiberguide.com

www.fiberguide.com




Figure 1: (a) Average SERS spectrum of seven fresh silver colloidal suspensions prepared with a batch process; (b) SERS spectrum of a potassium phosphate aqueous solution (final concentration of 5.0×10^{-3} M).

The objectives of this study were to unravel the characteristic SERS bands of *Salmonella typhimurium* bacteria and to verify the simple and universal algorithms established earlier (5). The ultimate goal is to develop label-free SERS for rapid, accurate, specific, and routine screening of *E. coli*, *L. monocytogenes*, and *S. typhimurium* bacteria for public safety and security.

Experimental

Chemical reagents. Chemical reagents (silver nitrate, >99%, and trisodium citrate, >99% from Sigma-Aldrich Co., St. Louis, Missouri) and tryptic soy broth (TSB) media (Becton, Dickinson and Co., Sparks, Maryland) were used without further purification.

Silver colloid preparation: During a 12-month period, 10 batches of citrate-reduced silver colloidal suspensions were prepared by a modified Lee and Meisel procedure (6) and were described in detail previously (5).

Bacterial culture preparation: More than 14 batches of *E. coli* strains (*E. coli* ATCC 25922), *L. monocytogenes* strains

(*L. monocytogenes* ATCC 13932), and *S. typhimurium* strains (*S. typhimurium* ATCC 14028) were incubated in TSB growth media at 37 °C for approximately 17–20 h without agitation in a 1-year span. This growth procedure routinely yielded a culture containing $\sim 10^9$ colony forming units (CFU)/mL of each respective bacterium at stationary phase.

SERS spectral collection and analysis. SERS spectra were collected on a Fourier transform (FT)-Raman module using a DTGS KBr detector and XT-KBr beam-splitter (Nicolet 670 FT-IR bench, Madison, Wisconsin). Raman scatter was accumulated using 180° reflective mode with 1 W of laser power at 1064 nm excitation and 256 scans at 8-cm⁻¹ resolution. All spectra were transformed into .spc files (Grams file format) and then were smoothed with the Savitzky-Golay function of two polynomial and 11 points by the use of Grams/32 software (Version 7.0, Galactic Industries Corp., Salem, New Hampshire). The data set was loaded into Microsoft Excel 2000 (Microsoft, Redmond, Washington) to execute simple algorithm analysis.




Figure 2: Representative SERS spectra of the same batch of *E. coli* in two batches of silver colloids. As a comparison, an FT-Raman spectrum of *E. coli* culture, labeled as intact, was provided. Spectra were shifted vertically for direct comparison.





Figure 3: Representative SERS spectra of the same batch of *L. monocytogenes* in two batches of silver colloids. As a comparison, an FT-Raman spectrum of *L. monocytogenes* culture, labeled as intact, was provided. Spectra were shifted vertically for direct comparison.

Two types of measurements were performed. First, 50 μL of different silver colloidal batches in both fresh and aging

state were introduced into 6 mm \times 50 mm commercial glass tubes (Fisher Scientific, Suwanee, Georgia) to evaluate



BALL MILLS

Retsch's ball mills are powerhouses with ultimate grinding efficiency and unparalleled versatility:

- 4 different planetary ball mills for ultra-fine grinding to sub-micron sizes
- 2 models of high-impact mixer mills for dry, wet, cryogenic grinding and cell disruption
- great variety of grinding containers, in sizes from 1.5 ml to 500 ml, in 7 different materials
- all-digital controls & memory for perfect reproducibility of results
- ease of use and outstanding safety features

GRAND PRIX 2008

Win with RETSCH at www.retsch.com/grandprix.

The first prize is a Mercedes sports car!



Retsch®

Solutions in Milling & Sieving

1-866-4-RETSCH

a VERDER company

www.retsch-us.com




Figure 4: Representative SERS spectra of the same batch of *S. typhimurium* in two batches of silver colloids. As a comparison, an FT-Raman spectrum of *S. typhimurium* culture, labeled as intact, was provided. Spectra were shifted vertically for direct comparison.




Figure 5: Plot of the ratio of I_{712}/I_{730} (horizon) vs. one of I_{390}/I_{352} (vertical) for *E. coli* cultures (\square), *L. monocytogenes* cultures (\bullet), and *S. typhimurium* cultures (\blacktriangle). Each point was the average of at least four SERS measurements in two silver colloidal batches that were stored for less than 40 days.

their reproducibility and stability, and 40 μL of an aqueous tribasic potassium phosphate solution was added to the same volume of silver colloids to examine the binding effectiveness of silver nanoparticles. Second, bacterial suspensions were mixed with silver colloids at a volume ratio of 1:1 (40 μL /40 μL) to get the respective SERS spectra. For one bacterial or potassium phosphate solution against one batch colloidal solution, two measurements were taken. Immediately after mixing two solutions, the tubes were shaken vigorously five times and then kept untouched for 10 min before the subsequent SERS measurements.

Results and Discussion

Assessment of reproducibility, stability, and binding effectiveness of silver colloids. One of the major concerns in SERS studies is to prepare reproducible silver colloids with the batch process, because the fabrication procedure is sensitive to many factors, such as temperature and rate of addition of reagents (7). As an alternative, the flow system has been developed to produce reproducible and stable colloids (8). However, batch-based colloidal fabrication is still attractive and practical because of its simplicity and its ease of operation (9).

Although UV-vis spectroscopy has been used to reveal the characteristics of silver colloids (3,7,8), Raman spectra can provide the information resulting from chemical residuals in colloids (5,9). Figure 1 shows the average spectrum of normalized SERS signals at 215 cm^{-1} from 10 fresh silver colloidal suspensions prepared with the batch process in the 1100–300 cm^{-1} region, in which no significant SERS peaks arising from decomposed chemical residuals in silver colloids were observed. Meanwhile, there were no additional SERS peaks observed for the silver colloids during an eight-week storage period at room temperature, suggesting the stability of the silver colloidal suspensions (5).

In addition, probe compounds such as dyes and inorganic phosphate have been introduced into colloidal solutions to assess their binding effectiveness (5,8). We selected potassium phosphate as an analyte because it is much smaller than

the silver colloidal nanoparticles and has characteristic P-O vibrations. Figure 1 (b) shows the SERS spectrum of an aqueous potassium phosphate solution at a final concentration of 5×10^{-3} M, and intense and separated ν_1 bands at 922 cm^{-1} were used to evaluate the binding effectiveness of silver colloidal nanoparticles at both fresh and aging states (5).

SERS characteristics of *E. coli*, *L. monocytogenes*, and *S. typhimurium* cultures: Figures 2 through 4 show typical SERS spectra in the region of $1100\text{--}300 \text{ cm}^{-1}$ from mixtures of two batches of silver colloidal suspensions and one batch of *E. coli*, *L. monocytogenes*, and *S. typhimurium* culture, respectively. As a comparison, FT-Raman spectra of individual cultures are given as intact, and they do not show any apparent bands. As expected, there are large variations in both relative intensity and position of SERS-active bands because of subtle changes among silver colloidal batches and unpredicted "hot" binding sites. Because incubated bacterial cultures consist of numerous species that were preexisting in growth media and were produced as by-products dur-

ing bacterial growth, clear understanding of the origins of SERS-active bands is not straightforward.

The most striking differences of SERS spectra among *E. coli*, *L. monocytogenes*, and *S. typhimurium* cultures are the appearance of more prominent SERS peaks in the *E. coli* and *S. typhimurium* suspensions than in *L. monocytogenes* suspensions. Such distinctions could be a reflection of their genetic, structural, or metabolic differences. Close examination of all SERS spectra over a period of more than one-year indicates that an intense 712-cm^{-1} SERS peak is common in *E. coli* and *S. typhimurium* suspensions and could be used to discriminate them from *L. monocytogenes* cultures, whereas the 390-cm^{-1} SERS band not only is indicative of *L. monocytogenes* cultures but also could be used to discriminate *E. coli* cultures from *S. typhimurium* cultures.

Identification of *E. coli*, *L. monocytogenes*, and *S. typhimurium* cultures: Ratios of I_{712}/I_{730} and I_{390}/I_{352} for three types of bacterial cultures were obtained by the same procedure we developed earlier (5) and are plotted in

Figure 5. Three cultures display much scatter distribution, probably as a result of a number of factors, such as the variations of chemical components and heterogeneities over batches of bacterial cultures and colloidal suspensions, and time after mixing bacteria and colloids. However, the greatness of the two ratios differs among the three species. Careful examination indicates that the *E. coli* and *S. typhimurium* cultures can be separated from *L. monocytogenes* by the ratio values of I_{712}/I_{730} greater than 1.04, and the *E. coli* and *S. typhimurium* cultures can be identified from each other by the ratio values of I_{390}/I_{352} around 0.97, respectively. Notably, the classification model is based upon SERS spectra acquired from different colloidal batches stored for different lengths of time.

Although multivariate data analysis of SERS spectra has been attempted to discriminate different bacterial species, SERS data from only one batch of colloid were used, and SERS spectra were preprocessed by subtracting a linearly increasing base-


μ MAX IR Microscope

- Sample Compartment Design
- Transmission, Reflection and ATR
- High Throughput Design
- Compatible with DLaTGS Detector
- Easy to Use

PIKE
TECHNOLOGIES

2901 Commerce Park Drive
Madison, WI 53719 USA
608.274.2721

www.piketech.com



Unravel the Matrix EDXRF Analyzers Jordan Valley AR Inc.

@Pittcon Booth 2134

Jordan Valley is a leading designer, manufacturer, and marketer of energy-dispersive X-ray fluorescence (EDXRF) systems and components with more than 30 years experience.

Our EDXRF instruments are unrivaled in quality, performance, flexibility and ease of use. Our powerful spectrometers are configured to specific industry requirements and designed to meet your targeted elemental analysis needs in the following markets:

PETROLEUM:

Ultra low-level sulfur detection

CHEMICAL:

Additive analysis & process control

PLASTICS/POLYMERS/WEEDS

RoHS/ELV:

Quantification of toxic metals

FORENSICS/HOME LAND SECURITY:

Mobile and on-site identification of unknown elements

ALLOYS:

Metal identification

MINERALS:

Major and trace element quantification

ENVIRONMENTAL APPLICATIONS:

Air particulate analysis and field screenings tool



8601 Cross Park Dr. Suite 200 Austin TX 78754
Toll Free: 1-877-777-1222

Visit our website at: www.jordanvalley-apd.com

line (3). Clearly, it is impractical and time-consuming. Consequently, development of ratio algorithms utilizing the unique SERS bands directly is simple and can be applied universally for fast, accurate, and routine screening of *E. coli*, *L. monocytogenes*, and *S. typhimurium* species on a variety of silver colloidal suspensions over batches and storages.

Conclusion

This study suggests the usefulness of the SERS technique in rapid and routine

species identification of *E. coli*, *L. monocytogenes*, and *S. typhimurium* cultures as a result of the technique's simple and easy sampling attribute as well as exclusive SERS biomarkers. To reveal consistent SERS bands from poor reproducibility of SERS signal of bacterial analyte, numerous spectra of mixing various batches of bacterial cultures with different colloidal batches were examined carefully. Characteristic bands at 712 and 390 cm^{-1} that appeared consistently and had the strongest intensity were identified. Two unique bands

were then used to identify *E. coli*, *L. monocytogenes*, and *S. typhimurium* cultures with 100% success.

To assess the reproducibility and stability of citrate-reduced silver colloids over batch process and over storage period as well as their binding effectiveness, FT-Raman spectra of silver colloid suspensions and ratio of SERS-active P-O band at 922 cm^{-1} in potassium phosphate aqueous solutions were used. The measures will ensure the quality control of batch-based silver colloids in the routine application over a long period.

Acknowledgment

Mention of a product or specific equipment does not constitute a guarantee or warranty by the U.S. Department of Agriculture and does not imply its approval to the exclusion of other products that may also be suitable.


References

- (1) C.W. Donnelly, in *Listeria, Listeriosis, and Food Safety*, 2nd ed., T. Ryser and E.H. Marth, Eds. (Marcel Dekker, New York, 1999), pp. 225–260.
- (2) T. Geng, J. Uknalis, S.-I. Tu, and A.K. Bhunia, *Sensors* **6**, 796–807 (2006).
- (3) R.M. Jarvis and R. Goodacre, *Anal. Chem.* **76**, 40–47 (2004).
- (4) M.L. Laucks, A. Sengupta, K. Junge, E.J. Davis, and B.D. Swanson, *Appl. Spectrosc.* **59**, 1222–1228 (2005).
- (5) Y. Liu, Y.-R. Chen, X. Nou, and K. Chao, *Appl. Spectrosc.* **61**, 824–831 (2007).
- (6) P.C. Lee and D. Meisel, *J. Phys. Chem.* **86**, 3391–3395 (1982).
- (7) C.H. Munro, W.E. Smith, M. Garner, J. Clarkson, and P.C. White, *Langmuir* **11**, 3712–3720 (1995).
- (8) R. Keir, D. Sadler, and W.E. Smith, *Appl. Spectrosc.* **56**, 551–559 (2002).
- (9) Y. Liu, Y.-R. Chen, K. Chao, and X. Nou, *J. Raman Spectrosc.* In press.

Yongliang Liu is with the Department of Nutrition and Food Science, University of Maryland, College Park, Maryland. **Yud-Ren Chen, Xiangwu Nou, Moon S. Kim, and Kuanglin Chao** are with the Food Safety Laboratory, Beltsville Agricultural Research Center, USDA, Beltsville, Maryland. ■

OIL-TRACE

Complete elemental analysis of fuels and oils



Oil-Trace overcomes the traditional challenges associated with XRF analysis of oils and fuels. XRF is acknowledged as a powerful technique throughout the petrochemical industry and Oil-Trace now takes it to the next level:

- dramatic reduction in number of calibrations
- significant operational cost saving
- automatic correction for dark-matrix composition
- analysis independent of variations in sample amount

Oil-Trace is a standard feature of Axios-Petro, and available as an upgrade to other PANalytical high-end WDXRF systems. The package includes all tools required for successful analysis according to international norms and methods.

For further information, contact your local representative or visit: www.panalytical.com

PANalytical Inc.
117 Flanders Road
Westborough, MA 01581
USA
T +1 508 647 1100
F +1 508 647 1115
Toll Free: +1 800 279 7297
info@panalytical.com
www.panalytical.com

