

A PCR- Based Assay Using Sequence Characterized DNA Markers for the Identification and Detection of *Aphanomyces euteiches*

Dr. George J. Vandemark et al.

USDA-ARS


Vegetable and Forage Crops Production Unit

Prosser, WA

Aphanomyces euteiches

Plant pathogen that causes severe root rot disease in alfalfa, peas, and beans

Affected Regions



Infected alfalfa



USDA-ARS, Prosser, WA

Limitations of Conventional Methods for the Detection of *Aphanomyces*

- Use of selective media is confounded by presence of other fungicide resistant microbes i.e. *Pythium*.
- “Baiting” soil with susceptible host requires up to three weeks for completion.
- Microscopic detection of oospores is tedious and oospores are only produced at the end of season.

Soil Baiting Technique with Aphanomyces



OBJECTIVES

- Design a system based on PCR that could discriminate *A. euteiches* from other closely related species and genera of soilborne microbes.
- Use the system to detect *A. euteiches* in infected roots.
- Use the system to detect *A. euteiches* in soil.

OBJECTIVES

- Design a system based on PCR that could discriminate *A. euteiches* from other closely related species and genera of soilborne fungi.
- Use the system to detect *A. euteiches* in infected roots.
- Use the system to detect *A. euteiches* in soil.

Advantages of PCR-Based Detection Assays for Soilborne Microbes

- Rapid
- Selectivity can be broad (genus) or narrow (species, race)
- Not necessary to isolate organisms in pure culture from soil/plant tissue
- Easier for ‘microscopically challenged’ pathologists

We chose to use SCARs for developing the assay

- SCARs = Sequenced Characterized Amplified Regions
- SCARs use a pair of sequenced characterized primers for amplification of target DNA
- SCARs were first developed for mapping resistance genes in lettuce to *Bremia lactucae*

Advantages of SCARs Over RAPDs


- Only detect a single locus (PCR product)
- Highly specific
- More rapid than RAPDs

Method for Designing SCAR Primers

- Identify a RAPD that is only amplified by all isolates of the target species (*Aphanomyces* sp.)
- Clone and sequence the RAPD
- Design primers based on RAPD sequence
- Optimize thermocycling reaction conditions.
 - [MgCl₂] (1.5 mM-4.5 mM)
 - annealing temperature (60°C-72°C)

SCARs

1 2 3 4 5



RAPDs

PRIMER gcctacactg



**CORTAR
FRAGMENTO**



**PURIFICACION
DEL FRAGMENTO
DE ADN**



**SECUENCIACION DE
EXTREMOS DEL
FRAGMENTO DE ADN**



**DISEÑO DE
PRIMERS
(20-25-mers)**

gcctacactgTTCCATGCATTACGG
gcctacactgGACGTAAGCTGATT





**AMPLIFICACION
DE ADN**

1 2 3 4 5



Developing a SCAR Specific for *A. euteiches*




BOTTOM

Effect of $[\text{MgCl}_2]$ on Amplification


1.5 mM MgCl_2

3.0 mM MgCl_2




Annealing Temperature = 60 °C

Optimized SCAR Reaction Conditions (1.5 mM MgCl_2 ; 70°C Annealing)





Effect of Cycle Number on Amplification¹

20 40 20 40




¹Two-Step PCR {94°C (1 min) ↔ 72°C (1 min)}

Developing a SCAR Specific for *A. cochlioides*



BOTTOM

A Single PCR Product (SCAR) is Diagnostic of *Aphanomyces euteiches*




Soil microbes tested with SCARs

<i>Aphanomyces euteiches</i>	21
<i>Aphanomyces cochlioides</i>	8
<i>Phytophthora infestans</i>	2
<i>Pythium ultimum</i>	3
<i>Pythium aphanidermatum</i>	1
<i>Pythium dissoticum</i>	1
<i>Fusarium oxysporum</i>	4
<i>Fusarium solani</i>	2
<i>Thelaviopsis basicola</i>	2
<i>Rhizoctonia solani</i>	2
<i>Mycosphaerella pinodes</i>	1
<i>Achlya</i> spp.	6


Multiplex PCR Demonstrates Species-specific Nature of SCAR Primers

Mixture *A. cochlioides* *A. euteiches*




SCAR Primers can Detect *A. euteiches* in Infected Roots


Infected roots H A. e



Detection of *A. euteiches* in Field Grown Plants



SCAR Primers can Detect *A. euteiches* in Organic Debris Fraction of Field Soil



Possible Applications for SCARs

● Qualitative

- Detection of pathogen in soil samples
- Detection of pathogen in infected tissue
- Discriminate between *Aphanomyces* spp.

● Quantitative

- Compare pathogen colonization/movement between different plant genotypes.
- Indirect selection for resistance among heterogeneous populations