



Glassware Washing Quiz

First Name:

Last Name:

Student/Staff ID (email):

Supervisor:

What should you do with chipped or cracked glassware:

Put in in the regular trash

Dispose it in the Broken Glass Waste Box

Dispose it in the Recycle Bin

Steam and heat are the most dangerous aspects of cleaning glassware.

True False

Filters and drains do not have to be checked first before loading glasses into the washing machine.

True False

Before loading containers of liquids into the autoclave, the caps must be loosened

True False

After autoclaving liquids, how long should you wait before you handle them?

10 minutes

30 minutes

45 minutes

1 hour 10 minutes

After autoclaving bottles, how long should you wait before you handle them?

5 minutes

10 minutes

30 minutes

Describe a good procedure to ensure uniform heating when autoclaving liquids

When acid washing, it is okay to keep bleach nearby.

True False

If you use chromic acid to clean glassware, you should

Ask your supervisor to select a safer cleaning method

Wear a heavy-duty apron and full-face shield when handling acid

Carry and store acid in a secondary acid-resistant container

Have another person stay nearby when you are using acid

NOTE! When finished, print out a copy, sign and date it, and return it to your Safety Officer

SIGN YOUR NAME HERE:

Date:

Last Modified: 04/04/2007