


NGS Applications for Diagnosis and Detection of Plant Pathogens




NGS Applications for Diagnosis and Detection of Plant Pathogens

- Brief Introduction to NGS Technology and Joint Projects
- NGS Capability at CFIA (Charlottetown Laboratory)
- Simultaneous Detection of Plant Pathogenic Bacterium and Viruses
- Summary and Future Prospects


Current Diagnostics Challenges in Plant Health



Brief History on DNA and DNA Sequencing



Big Frame on NGS Software and Databases



Computation Master Node Charlottetown Laboratory



Internal Storage System



Genome workstation



- Bacterial Refseq database
- Viral Refseq database
- Potato/Tomato Genomes
- Other reference databases
- Genome workstation • Genome workstation Genome workstation
- Genome workstation





Technolog Development

Diagnostics

Reference

NGS Hardware @ Charlottetown Lab



Ottawa Hospital Research Institute













NextSeq 500 @ AAFC Fredericton Lab

MiSeq @ Charlottetown Lab

Ion PGM @ AAFC Charottetown Lab



MiSeq & HiSeq

@ BC Cancer Reseach Center



MinION and GridION




PacBio System

NGS Software/Databases @ Charlottetown Lab



Window Based CLC Genomics Workbench



Computation Master
Node @ Charlottetown
Laboratory



Internal Storage System

Diagnostics

Genome workstation

Reference

- Bacterial Refseq database
- Viral Refseq database
- Potato/Tomato Genomes
- Other reference databases

Technology Development

- Genome workstation (Huimin Xu)
- Genome workstation (DL Hammill)*
- Genome workstation (Sean Li)
- Genome workstation (Donna Smith)*
- Genome workstation (Jingbai Nie)*
- Genome workstation (Eric Chuan)*

Linux-Based
Bioinformatics
Toolkits

AODP 2.4.6.4

Blast2Go

R 3.5.1

Mauve

CytoScape

BWA 0.7.17

BLAST+ 2.7.1

FastQC 0.11.8

AfterQC 0.9.7

SamTools 1.9

BBMap 38.26

BedTools 2.27.1

EDirect 9.90

SPAdes 3.12.0

GenomeTools 1.5.10

TrimGalore 0.4.5

Picard 2.18.14


FastXToolkit 0.0.14

AKATrim an in house

trimming tool

ANIGen Toolkit

Applications of NGS @ Charlottetown Lab



Research Group	Sample Type	Total No	Targets
	PPEQ	20	
Huimin Xu DL Hammill	DIA	10	Multiple Viruses with confirmation
	Research	19	
Sean Li JB Nie EJ Chuan	DIA	10	
	Imported	9	Bacteria and viruses with confirmation
	Research	15	
D. Smith	Potato wart	40	Transcriptomes
	GMO Trait	7	(wart) and artificial constructs (GMO)

RPS Project: CHA-P-1812A

Developing & synchronising methodologies with US toward high throughput detection and diagnostics of bacterial diseases

Project Lead: Dr. Sean Li, CFIA Charlottetown Lab, Charlotetown, Canada Collaborator: Dr. Michael Stulberg, USDA-APHIS, PPQ, Beltsville, USA

Objectives

- 1. To evaluate and validate the newly-developed and currently-used qPCR assays at CFIA Charlottetown Laboratory for detecting potato bacterial brown rot caused by *Ralstonia* solanacearum R3bv2 strains with methods employed by USDA-APHIS.
- 2. To carry out a comprehensive comparison of serological and molecular methods made available during the last 20 years for the detection and identification of *Clavibacter sepedonicus*.

CHA-P-1812A: Developing & Synchronising Methodologies with USDA-APHIS

Project Outcomes

- Generated reliable data for the accurate and sensitive detection of Rs R3bv2 and Cs in potato on the basis of currently available technology.
- Advise the Policy and Programs Branch of new developments in technology and implications for regulatory directives.
- Establish diagnostic schemes with automation component that will best fit in current diagnostic needs, and be transferred to the diagnostic labs of CFIA and other end users, e.g. private labs approved by CFIA.
- Enhance diagnostic capacity at CFIA and gain recognition of Canadian plant pathology expertise and research capacity.
- Publication in scientific journal for relevant research results and discoveries.


Other Contributors at CFIA with US Collaborators


- Plant Pathogenic Viruses: Drs. Mike Rott and Huimin Xu are the major contributors at CFIA. A software system, Virtool, has been established for practical applications under evaluation.
- Plant Pathogenic Nematodes: CFIA has bilateral connections and discussions with USDA-APHIS regularly.
- Potato wart: Dr. Andre Levesque and Donna Smith are the major contributors at CFIA, together with AAFC scientists.
- Dickeya-related potato blackleg disease: Dr. Sean Li is the major contributor in CFIA in collaboration with Drs. Amy Charkowski, Colorado State University, and Jianjun Hao, University of Maine.
- Bbr caused by Ralstonia solanacearum race 3 biovar 2 and Brr caused by Clavibacter sepedonicus: Dr. Sean Li is the major contributor in CFIA in collaboration with Drs. Michael Stulberg, USDA-APHIS, and Anne Alvarez, University of Hawaii.

Simultaneous Detection of PVS and CLso

- ➤ Potato zebra chip disease was named after its unique symptoms in potato tubers. It first appeared in Mexico in 1994, and was named in 2008 by Lia W. Liefting and her coworkers as 'Candidatus Liberibacter solanacearum'.
- The disease was initially observed as darkened bands in processing potato chips, due to an alteration in sugar metabolism.
 - ➤ Darkening of the medullary rays ranging from mild to severe in the entire length of the tuber.









Potato Zebra Chip Disease: the Pathogen and Vector

- 1. Potato zebra chip disease, also known as "papa manchada" or "papa rayada", is caused by the phloem-limited bacterium 'Candidatus Liberibacter solanacearum'
- 2. There are five haplotypes associated with the geographic ranges
 - Haplotypes A&B: Discovered in Mexico, North and Central US, and New Zealand, and vectored by Tomato/potato psyllid;
 - Haplotypes C,D &E: Discovered in Finland and North Europe, and vectored by carrot psyllid.
- 3. Vectors for the disease transmission include tomato and potato psyllid *Bactericera cockerelli* and carrot psyllid *Trioza apicalis*.



Potato/tomato Psyllids: Bactericera (Paratrioza) cockerelli




Carrot psyllid *Trioza apicalis*

Global Distribution of Potato Zebra Chip Disease



The Virulence of Potato ZC Pathogen



Genomics of Ca Liberibacter solanacearum

Tomato PCR-

Tomato PCR+


Potato PCR+

Potato PCR+ Potato PCR+

Potato PCR+

Potato PCR+ Potato PCR+





CFIA3-2 The tomato plant infected by haplotype A through grafting nine years ago.

This CLso haplotype A sample originated from potato plants from NDSU. Tomato plants and tissue culture plantlets were graft-transmitted from potato and preserved. All tomato plants graft-transmitted by haplotype B died.

WA5 The tomato plant infected by haplotype A through

The CLso haplotype A sample originated from potato plants from Washington. Tomato plants and tissue culture plantlets were graft

Genomic DNA Prep
















Genomics of Ca Liberibacter solanacearum



Simultaneous Detection of PVS and CLso

	CFIA3-2	WA5	C4	Carrot
'Ca Liberibacter solanacearum'	Haplotype A	Haplotype A	Haplotype D	Haplotype D
Unmapped Read Pairs	510043	1101999	235435	1870698
Mapped Read Pairs	5.72%	11.09%	3.71%	32.55%*
Potato virus S	Negative(PCR-)	Positive(PCR+)	Positive(PCR+)	Negative (PCR-)
Pepino virus	Negative(PCR-)	Negative(PCR-)	Negative(PCR-)	Negative(PCR-)

- Diagnostic conclusion: CLso haptotype A and D were identified in the four testing plant samples.
- Tomato plant WA5 originated from Washington State and C4 originated from Finland contain complete genomes of PVS, whereas tomato sample CFIA3-2 and Carrot sample contain no detectable level of PVS. RT-PCR confirmed the detection of PVS and CLso.
- All four samples contained small fragments homologous to Pepino virus (False Positive).
 PCR failed to detect any trace amount of the virus.
- Results were similar using both CLC Workflow and Linux Command Line analyses

Positive and False Positive

Tomato Plant WA5: Tomato plant grafted with potato with Clso haplotype A NCBI Blastn: E-value ≤ 1e-20; hit length ≥ 80 bp; Bit Score≥80)

Ca L. solanacearum	Potato Virus S	Ralstonia solanacearum
13 Contigs; BScore 179-717	17 Contigs; BScore 286-13K	22 Contigs; BScore ≤ 72
Hit Length: ≥ 218 (713) bp	Hit Length: ≥ 437 (7505) bp	Hit Length: ≤ 61 bp
Homology: ≥ 98%	Homology: ≥ 98%	Homology: ≥ 92%
Positive	Positive	False Positive

Tomato Plant C4: Tomato plant grafted with carrot with Clso haplotype D NCBI Blastn: E-value ≤ 1e-20; hit length ≥ 80 bp; Bit Score≥80)


Ca L. solanacearum	Potato Virus S	Ralstonia solanacearum	
3 Contigs; BScore 425-1379	6 Contigs; BScore 425-1399	13 Contigs; BScore ≤ 62	
Hit Length: ≥ 304 (430) bp	Hit Length: ≥443 (7505) bp	Hit Length: ≤ 61 bp	
Homology: ≥ 98%	Homology: ≥ 98%	Homology: ≥ 88%	
Positive	Positive	False Positive	

Metatranscriptomic Analysis (MiSeq vs NextSeq)

Features	CFIA3-2 Haplotype A		WA5 Haplotype A		C4 Haplotype D	
	NextSeq	MiSeq	NextSeq	MiSeq	NextSeq	MiSeq
Read Length	2x150	2x300	2x150	2x300	2x150	2x300
Total Reads	35.4M	17.9M	58.0M	20.1M	144.1M	12.7M
Mapped Reads	34.9M	11.1M	55.7M	17.9M	141.4M	12.4M
Unmapped Reads	526.9K	470.7K	2250.2K	2203.1K	2,752.6K	346.4K
% of Unmapped	1.49%	2.63%	3.88%	10.97%	1.91%	2.72%
Contigs/scaffolds	986	3,617	2344	15,238	1795	980
Blast hits (Clso)	6 (0.6%)	9 (2.5%)	5 (0.2%)	13 (0.09%)	3 (0.17%)	3 (0.3%)
Bit score (Clso)	376-1352	82-539	410-1386	179-717	372-470	425-1399

NextSeq: AAFC Federicton Lab; MiSeq: Charlottetown Lab

Analytic Outcomes of MiSeq Read Length



Final Identification: *Ca* L. solanacearum haplotype B



Metatranscriptome of B5 in MiSeq		Features	
2x150	2x300	Read Length	
10.6M	11.4M	Total Reads	
10.5M	11.1M	Mapped Reads	
128.3K	339.9K	Unmapped Reads	
1.20%	4.47%	% of Unmapped	
167	4,480	Contigs/scaffolds	
78 (46.7%)	386 (8.6%)	Blast hits (Clso)	
241-1206	212-1886	Bit score (Clso)	

Conclusional Remarks and Future Prospects

- Draft genome of Clso haplotype A was obtained through metagenomics approach. Using NGS and bioinformatics approaches, plants infected by bacterial and viral pathogens can be steadily detected simultaneously using NGS Techniques with qPCR confirmations.
- Tomato plant WA5 originated from Washington State and C4 originated from Finland contain complete genomes of PVS, whereas tomato sample CFIA3-2 and Carrot sample contain no detectable level of PVS. All four samples contained small fragments homologous to Pepino virus which were diagnosed as false positives. PCR assay failed to detect any trace amount of the virus.
- Results were confirmed using both CLC Workflow and Linux Command Line analyses Further analysis of differential trascriptomic expression of haplotypes A and B in tomato plants may demonstrate the key mechanisms in virulence and pathogenicity of haplotype B in tomato plants.
- The NGS and bioinformatics workflow are used in assisting diagnostic work on PPEQ and Diagnosis Section at Charlotettown

Acknowledgement

CFIA

Dr. Huimin Xu, Desmond L. Hammill, Jingbai Nie, Kat Yuan, Eric J. Chuan,
Donna Smith, Pamela Ross, and many others
Drs. André Lévesque, Mike Rott, Guillaume Bilodeau, Marc-Olivier Duceppe,
Burton Blais, Adam Koziol, and many others

AAFC Labs

Drs. André Lévesque, Wen Chen, James Tambong, Iyad Kandalaft, Xianzhou Nie, Kyle Gardner, Bourlaye Fofana

International Collaborators

Dr. Gary Secor @ NDSU, Dr. Joseph Munyaneza @ USDA-ARS, Dr. Mikko Lehtonen @ Finish Food Safety Authority Eviira,

Dr. Lia Liefting @ Diagnostic and Surveillance Services, Ministry for Primary Industries, New Zealand,

23 Drs. Charlie Rush, Fedeke Workneh, Li Paetzold @ TAMU

Thank You!